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X.C Manual for SpectraScope Spectral Data Acquisition & Processing Software for the S1000 Spectrometer

Save icons:

Click on these areas (top left of screen) to perform the same operations as:

- File -> Save -> Sample (rectangle icon filled with dots)
- File -> Save -> Dark (rectangle icon filled with black, make sure your setup has lamp off),
- File -> Save -> Reference (rectangle icon filled with white, make sure your setup has lamp on)

Data Cursor:

The data cursor is a vertical line which can be moved from left to right or right to left along the data scan. The data value at the cursor location is displayed in the top status bar.

Movement is controlled by selecting one of the four rectangles on top right of screen with arrows << < > >> using your mouse.

Alternately, the keypad keys will perform similarly:

Home (fast left) PgUp (fast right) End (left) PgDn (right) Del (remove cursor)

When using the mouse and pressing the right hand button, the data cursor will move to the next peak on the display, in the direction of the selected icon (left or right).

Panic icon:

Clicking this area (top of screen, just left of the data cursor icons) will cause the system to reset to average settings, in case you want instant relief from slow response.

It affects the ADC data sample frequency (-> 100 KHz) and the number of samples to average (-> 1).

System Indicators:

Observe the system indicators as yellow letters appearing in the lower left of the screen.

- y means that you have performed a y rescale
...to undo reselect View -> [desired mode]
- z means that you have performed an x zoom
...to undo reselect View -> x zoom off
- s means that you have performed a File -> Graph
...to undo reselect View -> Data in realtime
- r means that you have View -> Channel -> RatioMS
...to undo reselect View -> Channel -> Master

Interrupts:

The Adc (Analog to Digital Converter) adapter needs the use of an interrupt request (IRQ) in the PC plug in bus. SS defaults to IRQ7, however if you are having problems you may change this to IRQ3, IRQ4 or IRQ5 (try this first). At the Dos prompt type "SET SSIRQ=5" (if this works just put the same line in your AUTOEXEC.BAT file). If you have the PC1000 plugin spectrometer you must also set the switches on the board.

Speed up Episode:

The episodic data capture function initiated via the Setup -> Episode menu option writes continuous data to disk when there is no episode delay(0) or inter scan delays. Since your

disk has a rather large seek time as compared to rapid data capture requirements, you may need speed up this process. Install the Dos RAMDRIVE.sys command in your config.sys file.

```
Example: dos=high,umb
         device=\dos\himem.sys
         device=\dos\ramdrive.sys 2048 /e
         device=\dos\emm386.exe
```

NOTE: if you have MicroSoft Windows installed then use \windows\ in place of \dos\ above. This installs a 2 meg ramdisk in extended memory. This will hold approx 200 episode samples. You must tell SS where the ram drive is located by typing "SET SSRAMD=D" at the Dos prompt (also put this line in your AUTOEXEC.BAT file. "D" indicates that your ramdrive is drive "D". The ramdrive will always be the next letter above your last drive.

Fatal Error Codes	Error descriptions
01	Graphics driver register error
02	Graphics driver system error
03	Sample malloc (memory allocate)
04	Dark malloc
05	Ref malloc
06	Dbuf malloc
07	Sbuf malloc
08	Hbuf malloc
09	Write digitl(mode 13) [dout]
10	Config ctr0 (mode 10) [set_fdc]
11	Config ctr0 (mode 10) [set_fdc]
12	Load ctr1&2 (mode 17) [set_dsf]
13	Adc init (mode 0) [initDac]
14	Adc mux set (mode 1) [initDac]
15	DMA malloc [initDac]

File Menu Help summary:

New - allows the experiment name to be changed. The current experiment name is shown in the upper right status bar. The experiment name indicates the next filename to be used when performing a file save. The experiment name is the first six characters of the filename. The next two characters are the file number (00..99). The file number can be changed too. The file number is automatically incremented when doing successive file saves of the same file type. The file types cannot be changed as these are automatically selected for the context of the data. The file type is the last three characters after the period in the complete dos filename. Here are experiment files generated by SS when File->Save..options are used:

>>Format: data [, data] [if dual unit]
 eeeeeenn.cfg - contains setup configuration
 eeeeeenn.drk - dark data
 eeeeeenn.ref - reference data
 eeeeeenn.sam - sample data [when in scope mode]
 eeeeeenn.ep - episode data [via Setup->Episode]

File->Save->Sample generates a file type which depends on the current mode selected...

>>Format: wavelength, data [, wavelength, data]
 eeeeeenn.abs - absorbance mode data
 eeeeeenn.trm - transmission mode data
 eeeeeenn.rfs - ref spectra mode data
 eeeeeenn.ssm - sample data [scope mode]
 eeeeeenn.irr - sample data [irradiance mode]
 eeeeeenn.hc1 - history channel RED [any mode]
 eeeeeenn.hc2 - history channel GREEN [any mode]
 eeeeeenn.hc3 - history channel BLUE [any mode]
 eeeeeenn.hc4 - history channel YELLOW [any mode]

Open - allow the user to select various files to be read back in which were previously saved. These include an entire experiment (*.cfg, *.drk, *.ref, *.sam), or individually setup (*.cfg), or sample (*.sam), or dark (*.drk) or reference (*.ref). When combined with View -> Data graph STATIC selection, the *.sam files can be opened, graphed, and then plotted. Using this technique, sample data can be loaded, reprocessed (via setup options) and then saved as new samples in Scope, Absorbance and Transmission modes!

Save - allows the user to save data to files. This includes an entire experiment or selective files as described above in Open. Normally you would first Save -> Dark then Save -> Reference. These are required for the other view modes (ie: absorbance, transmission, ref spectra). The Dark and Reference data is not immediately saved to a file, so they can be repeated easily. This is different when selecting Save -> saMple, as it immediately generates *.ssm and *.sam in Scope mode! (* is the current experiment name and file number as indicated in the upper right status line). If you have selected Irradiance mode, then the file type generated is *.irr.

NOTE!!! Save -> saMple in other View modes (ie: absorbance, transmission, ref spectra), generate data files for that modes data type (*.abs,*.trm,*.rfs). These files can be looked at with File -> View. Optionally, these files can be graphed using File -> Graph menu, including samples (*.ssm) taken in Scope mode. The save icons on the top left of the screen perform the same functions as above, but very quickly: File -> Save -> (sample) (dark) (reference)

Plot - allows the user to plot the currently selected screen on the connected printer device. First use setup->printer configuration to select output resolution, output port, and type of print device. SS will remember your setup even after power off via the SS.cfg file. Plot provides an option to draw xy grid lines.

View - allows the user to look at the context of selected files.

Graph - allows the user to graph selected files on screen. These data files are:

*.abs (absorbance), *.trm (transmission), *.rfs (reference spectra), *.ssm (scope sample), *.ep (episode data) and *.hc* (history channel).

This selection may be successively repeated allowing many different files on the same graph. Only the first graph is automatically y_scaled. Each new graph overlay is selected as a different color.

When graphing episode data, many scans will overlay from a single file selection. Again each scan will become a new color. The colors will cycle over if the number of scan episodes exceed the number of available colors. If you select the 3D graph, each scan in the episode will be offset to fit the graph verticle axis. The scaling is set to accomodate all scans. When combined with the View -> X_zoom you can graph a selected range of the episode wavelengths. You must first select X_zoom the perform the File -> Graph -> Episode.

The New graph selection allows the next graph to first be cleared and then rescaled. This is automatic when repeating Episode graphs.

A fter a graph has been drawn, the system is left in a STATIC mode. To restore realtime operation, you must select View -> Data in REALTIME. A small "s" appears on bottom left of the graph to tell you the system is in STATIC mode. Also the mode name shown in the upper left part of the screen will indicate Mode: Graph.

Exit - back to where you came from (Dos or Windows).

Setup Help Menu Summary:

ADC sample freq: 1..500

Changes Analog to Digital Converter data sample frequency, in KHz. The ADC clock also runs the detector controller, and effectively sets the exposure or integration time for the CCD detector. The ADC rate should be adjusted for each experiment to maximize the detector output voltage and signal to noise ratio. The integration period is reported in a message box each time the ADC rate is changed. Please note that you should not attempt to set this value above the limit of your ADC device (125Khz for DAS16jr, 330KHz for DAS16-330, 500KHz for PC1000).

Range for ADC input: 1/2/4/8

Changes the input gain on the ADC device. It is normally set to 1=(0..10volts) to match the detector output voltage range, but can be changed to 2=(0..5v); 4=(0..2.5v); or 8=(0..1.25v) for amplification.

Number of samples: 1..9999

Sets the number of scans to signal average. Please note that the screen graph is updated AFTER this number of samples (scans) have been aquired.

Integration cycle: 1/4

Controls the CCD sensor integration period where 1=short (8192/2*ADC rate) or 4=long (32768/2*ADC rate). This is similar to the shutter control on a camera.

Optical pixel resolution: 1/5/9/17

This performs data smoothing by applying a moving average of adjacent pixels to the data arrays. For example, option 5 would average each pixel with the 2 pixels on the left & right.

Option 1 performs NO data smoothing. Smoothing levels 5, 9, and 17 correspond to the input fiber image size for 50,100 & 200 micron diameter fibers respectively.

Flash and scan timing:

S1000 spectrometers provide signals for triggering pulsed light sources such as xenon strobes at the D9 connector on the rear panel of the instruments. Pin D1 provides a middle of scan pulse only when data is being taken. Pin D6 provides a continuous train of pulses (1 or more per scan) regardless of data acquisition activity. These signals can be also used to synchronize (and trigger) external devices, such as lasers, with each spectrometer scan.

Inter scan delays can be setup to allow appropriate recharge time for flash units, or to adjust system timing patterns.

Correct for dark: ON/OFF

The first 16 "dead" pixels that are read from the CCD sensor give a useful approximation of the dark current even when the CCD is illuminated. If correct for dark is on, the mean values of the even and odd "dead" pixels are subtracted from the remaining even and odd pixels. This is typically accurate to within a count of "true" dark level and is especially usefull for comparing normalized data (see below).

Data normalization: ON/OFF

This mode divides observed voltages (in ADC counts) by the integration period (in seconds), allowing comparison of reference and sample data taken at different sampling rates.

AbsTrans drift correct: ON/OFF

Absorbance or Transmission mode drift correct option. This option can be used ONLY with a dual spectrometer set to the same wavelength range as the master. Changes in light levels detected by the slave are subtracted from the mastert, so as to compensate for drift in the light source.

View as differential:

Converts the spectral scans to the first derivative given as $d(\text{amplitude})/d(\text{wavelength})$. The result is a positive & negative trace about the zero axis. The derivative is approximated at pixel (n) as $(A(n+1)-A(n))/(wave(n+1)-wave(n))$.

Unit cal data & channel:

This allows the setting of the wavelength calibration coefficients for the spectrometer in use. In addition the user may select the physical Analog to Digital channel to attach the Master or Slave display to. This allows the Master and slave channels to be swapped. Pressing enter will make no changes. The coefficients are determined from a least squares fit to a second order polynomial. The intercept is the wavelength at detector pixel(0), the first slope is the dispersion in nm/pixels and the second corrects for the slight curvature. Spectrometers can be recalibrated with a wavelength standard such as our HG-1 low pressure mercury/argon lamp. All data arrays and files are based on detector pixels. The calibration coefficients are used to report the wavelengths associated with each pixel and for labeling graphs. The first, second, and third coefficients are taken from the calibration sheet supplied with each unit and may be input exactly as shown from the lines labeled in the following example:

SS prompt	Sheet label	example values
1st & 2nd coefs	X Coefficients(s)	0.577048 -5.792E-05
Enter Intercept:	Constant	226.025916

Episodic data sampling:

Performs unattended data save operations over a period of time. Each "episode" consists of signal averaging and processing the "number of scans" configured and by options in the SETUP and VIEW menus.

The number of episodes to be taken are entered along with the desired delay between episodes. The total time between edisodes is equal to the delay plus the time it takes to execute the episode. Since the mode (ie: scope, absorbance, transmission) controls the type of data taken (*.sam, *.abs, *.trm), you must remember to select the proper mode first. A "Suppress graph trace?" option allows the user to prevent display of each episode and speed up software execution.

History channel function:

Allows setup for the type of operations to be performed when the View->Channel->History is selected. There are upto 4 graph traces available: Red Green Blue Yellow. Each channel can be active or disabled. When active, the graph trace represents the resulting function value for each scan. The equation for the function is entered here, in the Setup->History channel function.

As the graph moves to the right, you may X_zoom & Y_rescale. Upon reaching the right hand graph border the entire graph will then be scrolled left as new data is added.

The history channel function library will save upto 10 different equations. Each of the 4 display channels can select equations from this library. A table of constants is provided for use in each equation. They are provided as symbols a1..a4 b1..b4 c1..c4 d1..d4.

Equations are entered using standard math notation for addition, subtraction, multiplication & division. Additionally: m(wavelength in nm) =master

s(wavelength in nm) =slave

b(pixel width) =band for wavelength

The above functions allow selection of several different wavelengths with different band widths, which could be used for specific process control applications.

eg: $b(5)*m(550)*a1+b(10)*m(600)*b1$

eg: $b(10)*(m(350)/s(350))*1.5$

In addition to standard math notations as shown above, any equation may use the following extended math functions:

ABS() ..absolute val
 ROUND() ..round up
 TRUNC() ..integer
 COS() ..cosine
 SIN() ..sine
 EXP() ..eponentiate
 SQR() ..square
 LOG() ..natural log
 SQRT() ..positive square root
 LOG10() ..log base 10
 TAN() ..tangent

Spectrometer type: single/dual

This should match the attached spectrometer hardware. However, if you wish to use just the master spectrometer in a dual unit, selecting single will speed up operations.

Screen Color Configuration:

Allows the default colors to be altered. A mouse is required to operate this panel!

Printer configuration:

Allows selecting the current output resolution for the File -> Plot function. The half low is the fastest while half med and half high are respectively slower. These are half page portrait plots. The land and full options are full page landscape and portrait plots in their respective resolutions low, med, and high. Some laser printers may not have enough memory to plot in full page resolutions.

You may also select the proper interface port and printer types. The following printers are supported:

HP LaserJet II and III, HP DeskJet bw and color,
Epson/IBM 9-pin dot matrix, Epson 24-pin dot matrix,
IBM Pro 24-pin, and IBM Quietwriter.

Note: most other brands of printers are compatible with one of these.

View Help Menu Summary:**Scope mode:**

Allows user to view "basic" spectrometer data in REALTIME or as a STATIC graph (see View->Data graph STATIC option below). Each pixel is graphed on the x_axis as wavelength in nanometers using the calibration parameters entered via setup. Reselecting View->Scope will clear any Y_rescale that had been previously selected. Scope data will be affected by options available in the Setup menu.

Absorbance:

Displays absorbance at pixel n using the current sample, reference, and dark data sets:

$$A(n) = -\log((\text{sample}(n) - \text{dark}(n)) / (\text{ref}(n) - \text{dark}(n)))$$

system indicator

Transmission:

Displays percent transmission at pixel n using the current sample, reference and dark data sets:

$$T(n) = ((\text{sample}(n) - \text{dark}(n)) / (\text{ref}(n) - \text{dark}(n))) * 100$$

Percent transmission is mathematically equivalent to percent reflectance and can be used when doing reflection experiments.

Irradiance:

Displays relative radiant energy at pixel n using the current sample, reference and dark data sets:

$$E(n) = (\text{sample}(n) - \text{dark}(n)) * B(n)$$

Where $B(n)$ is the computed component of the spectral distribution of the blackbody radiant emittance (at user selected temperature in degrees K) divided by the current reference(n) data set. The reference data set is established by observing a standard illuminant of known color temperature, and using the "Save Reference" function.

Ref spectra:

Displays the reference and dark data previously saved as ref(n)-dark(n). This is useful for detecting peaks in line spectra, etc.

Y_rescale:

When selected, Y_rescale evaluates the current data set for maximum values and then redraws the graph such that the current data will expand or contract to fit 80% of the Y scale. The Y scale is then relabeled to reflect the new magnitudes. This option can be successively repeated and is a significant tool when evaluating data in realtime. Please note the status indicator in the lower left corner of the screen "y". This indicates that Y has been rescaled.

X_zoom: ON/OFF

Allows a selected region on the X axis to be expanded to the full graph. This option will work on any type of graph, both STATIC and REALTIME, any selected mode (scope, transmission, absorbance, ref spectra), history channel and view differential. Please note the status indicator in the lower left corner of the screen "z". This reminds you that the Z axis has been zoomed. You must reselect View -> X_zoom to perform an "unZOOM".

Channel: Master/Slave/Both/History

Allows the user to view data from different spectrometers when using a dual unit. The "history" channel selection works for both single and dual spectrometers. It provides display of the history channel function which is configured via the Setup menu. For more detail, refer to the Help -> Setup -> History channel.

Qstatus: ON/OFF

The Quick status panel appears on the top right of the screen when turned on. It quickly summarizes the configuration selected in the Setup menu.

Graph trace: Dots/Lines/Overlay

Allows user to select graph types. When selecting overlay, the previously selected scan will NOT be erased. This provides interesting displays in REALTIME modes.

Data as: STATIC/REALTIME

REALTIME provides instant display of sample data processing configured via the Setup menu and View mode. When STATIC is selected, files may be reloaded and graphed. The file types are *.sam (samples taken in scope mode). Also *.ref (reference) and *.drk (dark) data may be loaded. These will affect the display modes for absorbance and transmission. The data can be manipulated via the Setup menu, reprocessed, displayed and saved again!

Please note the status indicator in the lower left corner of the screen "s". This reminds you that you are viewing STATIC data. After selecting File -> Graph the system is left in STATIC view. You must reselect this option to revert back to REALTIME viewing.

Apps Menu Help summary:**Auto calibration**

Calibration of spectrometer is performed automatically using a Mercury-Argon lamp when this option is enabled. This option is free with the purchase of the HG-1 calibration lamp.

Data communication

This option allows the spectrometer system to be controlled via the serial ports. Reflectance/Absorbance data can be exported to an external computer system on command.

Process ctrl Mgmt (2Q94)

The process control management option allows the operator to define upper and lower threshold values used to control discrete signals (0..5 volts) available on the PC printer port. These signals can be logically attached to the SpectraScope equation processor, which displays results in the red, green, blue, and yellow display channels.

Color Application (2Q94)

This option allows the user to derive standard L^* a^* b^* and xyz tri stimulus values for color samples in realtime graphs in the SpectraScope environment.